

# MOISTURE ANALYZER

## ALH MB 1



# Moisture Analyzer

ALH MB 1




## MOISTURE ANALYZER

ALH MB 1

### SALIENT FEATURES

- Halogen Light Heating
- Touch Screen, Loadcell Sensor
- Temperature and time can be set
- The percentage of Moisture content
- Curve (Graph) residual
- The percentage of dry
- Store historical setting
- RS232C Interface



#### MOISTURE DETERMINATION (1)

Machine Type : PGB1MB  
Heating Mode : STANDARD  
Stop Mode : AUTO STOP  
Heating Temp. : 130°C  
Elasped Time : 03:25  
MW : 2.056g  
DW : 1.933g  
Moisture : 5.98%M  
Date : 2017-10-17  
Time : 08:16:14



**STANDARD ACCESSORIES :** Main Instrument, Aluminium PAN (50nos), Dust Plate, PAN Handle, PAN Supporter, Power Cable, Instruction Manual, Test Certificate.

## TECHNICAL SPECIFICATION

HSN Code - 90160020

MODEL	ALH MB 1
Make	Althea
Capacity	110g
Readability	0.001g ( 1mg )
Repeatability ( ± )	0.001g
Linearity ( ± )	0.002g
Display	Touch Screen
Pan Size	90 mm Ø
Sensor	loadcell
Calibration	Automatic External
Operating Temperature	5°C to 35°C
Heating Source	Halogen lamp
Temp. Sensor	PT - 100
Interface	RS232C
Temperature Setting	40°C ~ 199°C
Temp. Step	1°C
Time Setting	1~99 minutes
Dry Readability	0.01%
Dry Range	100.00% ~ 0.00%
Curve	Yes
Aluminium Sample PAN	Available
No. of Storage	20
Power Supply	AC 220V / 50-60 Hz
Dimension (LxWxH)	490 x 360 x 350 mm
Weight	7.5 kg

Specifications are subject to change without notice.