

All-metal body to create

Use high-quality metal material,
hard quality is the key.



- Food & Beverage Industry
- Pharmaceutical & Life Sciences
- Cosmetic Industry
- Environmental Testing
- Clinical & Diagnostic Labs
- Toxin Analysis

Safety Door

Mirror stainless steel
Special toughened glass
Beautiful and not easy to deformation,
greatly enhance the ability of anti-collision.

STOMACHER BLENDER

3/12 TIMES PER SECOND, HIGH EFFECTIVE WORKING

FOOD




COSMETICS



MEDICINE



TOXIN ANALYSIS




**HIGH QUALITY METAL
ELEGANT DESIGN**



**MIRROR
STAINLESS
STEEL**

+

Video dictionary
desktop support fu
l-text translation, smart
adjustive
Nt order, support
screen words, quick
click download



**EASY OPERATION
BLENDER STOMACHER**

MICRO COMPUTER MULTISTAGE PROGRAMMING



Features

Althea allows sterile homogeneous solid samples extracted from the bacteria process becomes very simple, just to sample and diluent to sterile sample bags, sample bags and then placed in the slap-type homogenizer to complete the processing of samples. Effective separation of the solid samples were included in the internal and surface micro-organisms homogeneous sample to ensure that all the samples sterile mixed bag. Treated sample solution can be sampled and analyzed, a sample of the change and the risk of cross contamination.

Key feature

- Large-screen LCD
- Adjustable homogeneous time.
- Homogenization speed adjustable or fixed.
- Slap before and after the device adjustable distance.
- Sterile disposable bag, to ensure health and safety.
- All open doors for easy cleaning
- Easy to observe the transparent window gla
- Samples with homogeneous non-contact instrument, if no sample leakage is not required for system cleaning.
- Homogeneous soft, sample no pollution, no damage, no heating, no sterilization, without washing.

Disinfection, Built-Indisinfection Lamp

Homogenization process can disinfection samples and improves work efficiency.



O3
Sterilization
without dead ends



Ultraviolet
sterilization



Flapping



Safety
monitoring



NTC
Detecting
system



Parameters are as below

Model	ALHSB-400	ALHSB-1000	ALHSB-3750
Mixer Type	Paddle		
Volume	3-400ML	50-1000 ML	200-3750 ML
Paddle speed range	180-720 RPM	360-1440 RPM	
Paddle speed	3~12 times of squeezing/second		
Paddle spacing	0~55 mm (Adjustable)		
Timer Range	1s-99h59m59s (8 Set of timing values can be programed simultaneously)		
Control mode	Microcomputer control		
Parameter storage	3 segment combination		
Disinfection	Yes with UV Light function		
Display	LCD Touch Screen		
Display Mode	4.3 inch touch screen display		
Homogenizer Bag	50 Bags	10 Bags	10 Bags
Overall Dimension	420*230*330mm	420*430*440mm	425*450*440mm
Machine weight	20 KG	36 KG	42 KG
Sterile bag size	170*300 mm	300*500 mm	355*500 mm
Cabinet material	Stainless steel +anti-corrosion coating		
Start mode	soft start		
Pause function	YES		
Power supply	Power : 220V/50hz-220w		
Anti-pinch function	Automatic stop,anti-pinch		
Other functions	Load and unload window + toughened glass door		
Rattling speed	3~12 times/sec (to 1 times/second step) Display the actual rattling number, the system automatically adjusts the slap speed Multi-Stage programming can be designed: parametric groups can be set according to one's wishes.		
End Alarm function	Yes		
Sterilizing Function	Yes		

DETAILS SHOW THE QUALITY

SURFACE PROTRUDING

The Protruding Starting Point On The Striker Plate Makes The Homogenization Better. Thickening And Longer Service Life



HEAT DISSIPATING HOLE

The radiator installed on the 3 sides of the fuselage to ensure that the machine does not heat easily for a long time



DETACHABLE DOOR HANDLE

304 stainless steel is built with longer service life.

The homogenizing effect is better for the homogeneity of the sample



STERILE BAGS 50 PCS

The homogenizing effect is better for the homogeneity of the sample



FUSELAGE IS MADE OF ALL METAL

52 LANGUAGE TRANSLATION, SUPPORT VOICE TRANSLATION AND PHOTO TRANSLATION

AESTHETICS

ANTIOXIDANT

LONG LIFE

STABILISATION

Warranty

- 12 months warranty for Althea Stomacher Blender.
- If any technical problem during warranty time, replacement parts will be sent free of charge.
- Technical support is also available in whole time.

Our Advantage

- Reply your enquiry in 24 working hours.
- Experienced staffs answer all your questions in fluent English.
- OEM is available if the quantity is large.
- We offer very competitive for all our products.

After-sales

Our machines undergo an extensive and strict testing and quality control procedure, ensuring every item purchased meet global quality standards. It's safe and reliable quality for Stomacher Blender.

