

MICROPLATE ELISA WASHER

ALH EW-1261



MICROPLATE ELISA WASHER

MODEL : ALH EW-1261

PRODUCT FEATURES :

- Large LCD display, intelligent software, easy to operate
- Strong washing function, bottom washing and two positions aspiration make completed washing effect
- 8-channel or 12 channel optional
- Needle positions adjustable and up to 10 different needle positions pre-settable Pressure, volume & time adjustable, washing times, soap time pre-settable Pressure, volume & time adjustable, washing times, soap time pre-settable, whole plate or single strip reading
- Self calibration and trouble alarm function
- Large memory to store up to 100 program



ACCESSORIES :

- 4 wash bottle (2500ml)
- Connecting tube set
- Cleaning wire

- Main lead
- Instruction Manual



MICROPLATE ELISA WASHER

MODEL : ALH EW-1261



TECHNICAL SPECIFICATION :

Plate	48/96 well and trips (U.V. Flat bottom)
Residual Volume	≤2µl
Dispensing volume	Adjustable from 50 to 900 µl
Wash cycles	1-16 optional
Mode of wash	From 1 to 12, row and rank, selectable
Soak time	Up to 990 seconds adjustable
Shake time	Up to 990 seconds adjustable
Display	Large LCD
Channel	3 channels for washing liquid, 1 channel for waste liquid
Memory	up to 100 programs can be stored
Interface	RS232
Dimension	450(L)mm x 310(W)mm x 170 (H)mm
Net Weight	7.5kg
Power	AC220V±10%, 50/60Hz; AC110V±10%, 50/60Hz



Althea Solutions Pvt. Ltd. GST PROVISIONAL ID :- 27ABDCA9914G1ZK

Address : 142, First floor, Lodha Boulevard, Majiwada, Thane (W)-400601, Mumbai, India

Factory : R-68 T.T.C. MIDC Industrial Area, Rabale, Navi Mumbai, MH- 400701, India

Mobile : +91- 9930833683 / 8928808788

Email : sales@altheatechnology.in/ info@altheatechnology.in/ milan@altheatechnology.in

Web : www.altheasolutions.in